

# Pharmachine INDIA

## **Co-Mill**



More Accurate Particle Size Control with Increased Output



## Co-Mill


Co Mill is used for Dry and Wet milling/size reduction of wide range of products. It is suitable for stand alone as well as online milling.

#### **Benefits & Highlights**

Low speed design allowing for inherently gentle action on the product resulting in optimum particle granularity, low heat, low noise, minimal fines, consistent, high quality end product, can be used on heat sensitive products, high capacities, fully customizable to many different applications and process stream.

### Advantages as compared to other Milling **Equipment:**

- The Output in Co Mill is very near to the size of the sieve used approx 80% - 90% where as in Multi Mill it is 1/3 of the sieve
- The Fines generated are on a lower side as compared to other mills
- · Heat generation in Co Mill is comparatively very low as it works on forced pressure
- · Other Mills have more of dusting as compared to Co Mill
- · Output in Co Mill is more than any other
- Design is flexible as per your requirement.







FAX: +91(22) 2577 1900

E-MAIL: pi@pharmachine.com

WEBSITE: www.pharmachine.com

52, SUYOG INDUSTRIAL ESTATE, LBS MARG. VIKHROLI (W), MUMBAI - 400 083, INDIA